

Breathable
Heavy-or-
special-
patients



The Turning Sheet is used for pain-sensitive patients, heavy patients. Moving and laterally repositioning patients is made much easier. Rolling up the Turning Sheet against the patients body in a lateral position and then fixing it in place with a bed sheet provides simple and gentle lateral position support.

The Turn Sheet is very versatile. Depending on the type of restriction, weight and sensitivity to pain we recommend using it as a complete set PM-1580

Technical specification

Product Name	Material	Colour	Max Lifting weight	Guarantee	Storage/transport
Turn Sheet	Cotton, natural thickening, latex free	Cream	200 kg	3 years	5 °C – 40 °C protect against sunlight

Item no.	Dimensions in cm	Weight	UDI
PM-1582	140 x 200	1950 g	(01)04250355315823(10)LOT

Hygiene and maintenance

The product doesn't require any special maintenance. It can be machine washed at 90°C and dried at low temperatures (45° – 50°C). We recommend washing with mild detergents. The use of cleaning agents with bleaching ingredients or other chemical laundry additives, also as fabric softeners may lead to permanent deterioration of the product properties.

